VIRUSOLVE®+ SPORICIDAL WIPES

SPORICIDAL DISINFECTION & CLEANING FOR MEDICAL DEVICES AND ANY SURFACE

Virusolve®+ is a patented proven HLD (High Level Disinfectant). The formulation does not contain any hazardous aldehydes or chlorine generating components.

Virusolve®+ is highly effective as a cleaner and terminal HLD against bacterial spores, mycobacteria, bacteria, fungi and viruses. The product is designed to eliminate the risk of cross-infection from spores, TB, viruses, fungi and bacteria.

Virusolve®+ does not simply render the microorganisms inert but has been demonstrated to kill these by disrupting the RNA of the microorganisms. It is safe to use on any surface and is odourless and colourless.



Primary features attributable to Virusolve®+ are:

- High Level Disinfection and cleaning.
- Does not contain chlorine or aldehyde compounds.
- Effective against bacterial spores e.g. Clostridium difficile (C dif), Bacillus cereus etc.
- Effective over a wide range of pathogens including mycobacteria, viruses, bacteria, fungi, yeasts and moulds, including MRSA, Hepatitis B and HIV.
- Effective against the difficult non-enveloped type viruses such as Polio, Adenovirus and Norovirus.
- Non-fixative.
- Non-corrosive will not damage delicate coatings.
- · No activation required.
- Cost effective.
- No fumes.
- Tested to EN14885 for disinfection properties

Packaging

- Ready to Use
- 200 soft pack
- 225 tub
- 225 bucket

Guaranteed protection and performance, giving total infection control - safely



DISINFECTANTS

Application

Virusolve®+ sporicidal wipes are ideal for removing biohazards in hospitals, health centres, dentists, vets, nursing homes, laboratories, food preparation areas and any environment where the risk of cross-infection needs to be eliminated, e.g. cruise ships, disinfecting hairdressing equipment etc.

Virusolve®+ sporicidal wipes are ideal for the cleaning and disinfection of hard surfaces, walls, floors and glasswork.

Virusolve®+ sporicidal wipes are ideal for HLD (High Level Disinfection) of medical devices and surgical instrumentation.

Directions for Use

Pre-impregnated Virusolve®+ sporicidal wipes are supplied in containers and are ready to use direct from the container.

General Surface Cleaning and Healthcare Environment:
Take a wipe and wipe gently over the whole area to be cleaned. Allow a minimum contact time of 1 minute before wiping clean to remove grease, blood and organic debris.
Removal of stubborn deposits may be assisted by the use of a stiff bristle brush. Rinse treated surface with clean water if required. Vinyl surfaces must be rinsed to avoid discolouration.

Note: Food preparation surfaces should be allowed to dry before use.

For Disinfection of medical, surgical and dental instruments and equipment (Manual procedure):

Wipe all surfaces of the instrument (preferably disassembled, where applicable) using a pre-impregnated Virusolve+ Wipe and allow the fluid to remain in contact for at least 5 minutes.

Examine the instrument and wipe or brush clean to remove debris that remains.

Rinse treated instrument in sterile, cold water and allow to dry.

Replace used Virusolve®+ wipes frequently or when gross soiling of the wipe becomes evident.

VIRUSOLVE®+ SPORICIDAL WIPES SPORICIDAL DISINFECTION & CLEANING FOR MEDICAL DEVICES AND ANY SURFACE

Virusolve®+ Properties

Active Against	Standards	Contact Time
Spores	EN13704 EN14347	2 Minutes 1 Minute
Mycobacteria	EN14348 Mycobacterium terrae Mycobacterium tuberculosis EN14563 Mycobacterium terrae Mycobacterium avium	1 Minute 1 Minute 5 Minutes 5 Minutes
Bacteria	EN1276, EN13727, EN13697, EN14561, EN16615	1 Minute 5 Minutes
Viruses	EN 14476 Polio, Adeno & Noro HIV, Hep B & C	5 Minutes 1 Minute
Fungi / Yeasts	EN1275, EN1650 EN13624, EN14562, EN16615	1 Minute 5 Minutes





Amity Limited, Libra House, West Street, Worsbrough Dale, Barnsley, South Yorkshire, S70 5PG, UK

T: +44 (0) 1226 747660 **F:** +44 (0) 1226 291936 **E:** sales@amityinternational.com

WWW.AIHC.CO.UK

© Amity Limited all rights reserved 2018

DISTRI	BUTE	D BY
--------	------	------