



**MODELLO: SONOACE 8000 Live**

**PRODUTTORE: MEDISON**

### **Specifiche tecniche**

- Color Doppler Imaging and Angio
- Pulsed Wave (PW) Spectral Doppler
- Tissue Doppler Imaging
- Tissue Harmonic Imaging
- Pulsed Inversion Harmonic Imaging
- Full digital beamforming
- Multi-frequency steps (General, Resolution, Penetration)
- 256 shades of gray
- Cine for 128 frames and Loop Review for 2048 lines
- Color VGA monitor with non-interlaced display
- OSIO(Auto image optimization)
- Auto Measurement
- SonoView II Image management
- Combined modes:2D/M, 2D/PWD, 2D/C, 2D/CPA, 2D/Steered CW, 2D/CPA/PWD, 2D/C/PW(TDI) , Steered CW, 2D/C/PWD Integrated 3D Imaging Package
- VOCAL(3D volume measurement package)
- MagiCut
- Weight: 222 lbs. / 101kg
- Height: 59" / 149 cm
- Depth: 34" / 85 cm
- Width: 26" / 67 cm