




SAMSUNG SMART
DIGITAL RADIOGRAPHY



XGEO GR40W

INNOVATIVE WIRELESS DETECTOR
FOR YOUR DIAGNOSTIC CONFIDENCE



INNOVATIVE WIRELESS DETECTOR, XGEO GR40W OFFERS CONVERSION OF ANALOGUE/CR INTO DIGITAL RADIOGRAPHY WITH EXCEPTIONAL CONVENIENCE

Introducing the new "Innovative Wireless Detector, XGEO GR40W" by Samsung Electronics! Features include wireless convenience, light and compact design, and high quality image preview (within 3 seconds).

The XGEO GR40W features Samsung Electronics state-of-the-art, advanced technology, which provides strong diagnostic confidence at a competitive price and allows for easy upgradability from analog X-ray or CR to DR.

Basic Units of XGEO GR40W

35 x 43cm Csl Portable

XGEO Station



Samsung XGEO GR 40W

User Friendliness



Compact Design & Rapid Use



Premium Image Quality at High Speed



High Connectivity



Samsung MoVue™



EXCEPTIONAL DIGITAL RADIOGRAPHY EXPERIENCE WITH XGEO GR40W

Portable, Easy, & Versatile Wireless Detector

The XGEO GR40W is the same size as a film, CR cassette making it fast and easy to convert analog X-ray or CR into DR. Innovative Wireless technology provides the ultimate convenience.

Compact Design, Simplified Workflow, & High Efficiency

Customer satisfaction will be high due to the compact design (35X43cm & 2.9kg) and simplified workflow. Efficiency is improved because the battery time is up to 30 minutes in battery charger mode, and up to 60 minutes in tethered mode, allowing the detector to be operated continuously.

Exceptional Image Quality & only 3 Seconds Preview Time

The XGEO GR40W provides exceptional image quality while minimizing dose to the patient. Industry leading post processing technology allows for image preview in only 3 seconds!



High Connectivity

Quick diagnosis improves user throughput, and ensures the best communication with patients.

Samsung MoVue™ *

We provide the portable Samsung MoVue™ to carry images and interact with patients more efficiently. It is applicable to medical imaging such as Ultrasound (3D), CT, and MRI. Samsung MoVue software is provided for Galaxy Tab to maximize diagnosis efficiency.

+ PACS + DICOM 3.0 compatible



* Option